


# PerioType<sup>®</sup> X-Pert

## You can clean it like a tooth!

*In the past abutments scratched easily during regular cleaning, causing rapid adhesion of bacteria.*

*The new patented zirconium nitride ceramic coating of the abutment makes cleaning for the dental hygienist easy. They may use the same instruments, and treat the implants just like real teeth.*

*Treating your patients with zirconium nitride coated abutments can give you the confidence for long-term success.*



### ZircoSeal<sup>®</sup>, zirconium nitride ceramic coating for abutments

- ✓ Is 3 times harder than steel currettes used to clean teeth and implants
- ✓ Is 6 times harder than titanium
- ✓ Has proven in vitro to produce less plaque adhesion
- ✓ Has proven in vitro to produce less pathogenic bacteria
- ✓ Has been proven to cause less inflammation than titanium implants



Titanium abutment  
after use of currettes

ZircoSeal<sup>®</sup> stays  
untouched by currettes